



Mouse anti-Human JRKL monoclonal antibody, clone 5F21 (CABT-B10485)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	JRKL (NP_003763, 426 a.a. ~ 525 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	5F21
Conjugate	Unconjugated
Applications	WB, ELISA
Sequence Similarities	ENLEKWLEVDSTEPGYEVLTDSEIIRRAQQGADESSSENEEEEEIELIPEKHINHAAALQWT ENLLDYLEQQGDMILPDRLVIRKLIRATIRNKQKMTKSSQ*
Format	Liquid
Size	50 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The function of this gene has not yet been defined, however, the encoded protein shares similarity with the human (41% identical) and mouse (34% identical) jerky gene products. This
---------------------	---

protein may act as a nuclear regulatory protein. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, May 2012]

Keywords	JRKL; JRK-like; HHMJG; jerky protein homolog-like; human homolog of mouse jerky gene protein;
-----------------	---

GENE INFORMATION

Entrez Gene ID	8690
-----------------------	----------------------

UniProt ID	Q9Y4A0
-------------------	------------------------

Function	DNA binding
-----------------	-------------
