



Mouse anti-Human KBTBD10 monoclonal antibody, clone 2I4 (CABT-B10490)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KBTBD10 (NP_006054, 205 a.a. ~ 305 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	2I4
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	RTDKENRVKNLSEVFDCIRFRLMTEKYFKDHVEKDDIIKSNPDLQKKIKVLKDAFAGKLP EPSKNAAKTGAGEVNGDVGDEDLLPGYLNDIPRHGMFVKD*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Kelch repeat and BTB domain-containing protein 10 is a protein that in humans is encoded by the KBTBD10 gene. Mouse monoclonal antibody raised against a partial recombinant
---------------------	--

KBTBD10.

Keywords

KBTBD10; kelch repeat and BTB (POZ) domain containing 10;

GENE INFORMATION

Entrez Gene ID

[10324](#)

UniProt ID

[O60662](#)
