



# Mouse anti-Human KIAA1967 monoclonal antibody, clone 5F7 (CABT-B10516)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                                                                                                                                                                   |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>             | KIAA1967 (AAH18269, 5 a.a. ~ 366 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.                                                                                                                                                            |
| <b>Isotype</b>               | IgG1                                                                                                                                                                                                                                                                              |
| <b>Source/Host</b>           | Mouse                                                                                                                                                                                                                                                                             |
| <b>Species Reactivity</b>    | Human                                                                                                                                                                                                                                                                             |
| <b>Clone</b>                 | 5F7                                                                                                                                                                                                                                                                               |
| <b>Conjugate</b>             | Unconjugated                                                                                                                                                                                                                                                                      |
| <b>Applications</b>          | WB,IF,sELISA,ELISA                                                                                                                                                                                                                                                                |
| <b>Sequence Similarities</b> | LPEKVVSPPEPEKEEEAAKEEATKEEEAIKEEVVKEPKDEAQNEGPATESEAPLKEDGLLP<br>KPLSSGGEEEEKPRGEASEDLCEMALDPELLLLRDDGEEEFAGAKLEDSEVRSVASNQSE<br>MEFSSLQDMPKELDPSAVLPLDCLLAFVFFDANWCGYLHRRDLERILLTLGIRLSAEQAK<br>QLVSRVVTQNICQYRSLQYSRQEGLDGGLPEEVLFGNLDLLPPP GKSTKPGAAPTEHKAL<br>VSHNGSLINVGSLLQ |
| <b>Format</b>                | Liquid                                                                                                                                                                                                                                                                            |
| <b>Size</b>                  | 100 µg                                                                                                                                                                                                                                                                            |
| <b>Buffer</b>                | In 1x PBS, pH 7.2                                                                                                                                                                                                                                                                 |
| <b>Storage</b>               | Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.                                                                                                                                                                                                          |

## BACKGROUND

|                 |                                                                                                                                                                          |
|-----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Keywords</b> | DBC.1; DBC1; deleted in breast cancer 1; DBC-1; p30 DBC protein; NET35; p30DBC; Deleted in breast cancer gene 1 protein; p30 DBC; OTTHUMP00000225106; OTTHUMP00000225109 |
|-----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

---

## GENE INFORMATION

|                       |                       |
|-----------------------|-----------------------|
| <b>Entrez Gene ID</b> | <a href="#">57805</a> |
|-----------------------|-----------------------|

---

|                   |                        |
|-------------------|------------------------|
| <b>UniProt ID</b> | <a href="#">Q8N163</a> |
|-------------------|------------------------|

---

|                |                                                                |
|----------------|----------------------------------------------------------------|
| <b>Pathway</b> | TNF-alpha/NF-kB Signaling Pathway, organism-specific biosystem |
|----------------|----------------------------------------------------------------|

---

|                 |                                                                                       |
|-----------------|---------------------------------------------------------------------------------------|
| <b>Function</b> | enzyme binding; enzyme inhibitor activity; enzyme inhibitor activity; protein binding |
|-----------------|---------------------------------------------------------------------------------------|

---