



Mouse anti-Human KIF26A monoclonal antibody, clone 4D21 (CABT-B10519)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KIF26A (AAH09415, 1 a.a. ~ 204 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	4D21
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MSESGAASPGARTRSLKSPKKRATGLQRRRLIPAPLPDTTALGRKPSLPGQWVDLPPPLA GSLKEPFEIKVYEIDDVERLQRPRTPREAPTQGLACVSTRLRLAERRQQRLREVQAKHK HLCEELAETQGRLMLEPGRWLEQFEVDPELEPESAEYLAALERATAALEQCVNLCKAHVM MVTCFDISVAASAAIPGPQEVDV*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords KIF26A; kinesin family member 26A; kinesin-like protein KIF26A; KIF26A variant protein;

GENE INFORMATION

Entrez Gene ID [26153](#)

UniProt ID [Q9ULI4](#)

Pathway Factors involved in megakaryocyte development and platelet production, organism-specific biosystem; Hemostasis, organism-specific biosystem; Kinesins, organism-specific biosystem

Function ATP binding; NOT ATPase activity; microtubule binding; NOT microtubule motor activity; nucleotide binding; protein binding
