



Mouse anti-Human KIR3DX1 monoclonal antibody, clone 4F8 (CABT-B10524)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KIR3DX1 (AAH33195.1, 1 a.a. ~ 78 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	4F8
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MANTEPTGQRTDEEEPAAEETQEIIYAQLNHQALSQTGFPPASQCPhyLSEDPSIYITV HQAQAEARAAPSLWHKGH
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords KIR3DX1; killer cell immunoglobulin-like receptor, three domains, X1; LENG12; KIR3DL0;

GENE INFORMATION

Entrez Gene ID

[90011](#)
