



Mouse anti-Human KIRREL2 monoclonal antibody, clone 3C02E4 (CABT-B10525)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KIRREL2 (AAH07312, 1 a.a. ~ 232 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3C02E4
Conjugate	Unconjugated
Applications	WB,IP,sELISA,ELISA
Sequence Similarities	MLRMRVPALLVLLFCFRGRAGPSPHLLQQPEDLVLLGEGGAQASLGRRASASFSEQKNL MRIPGSSDGSSSRGPEEEEETGSREDRGPIVHTDHSDLVLEEEGTLETKDPTNGYYKVRGV SVSLSLGEAPGGGLFLPPPSPLGPPGTPTFYDFNPHLGMVPPCRLYRARAGYLTPHPRA FTSYIKPTSFGPPDLAPGTPPFPYAAFPTPSHPRLQTHV*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Kin of IRRE-like protein 2 is a protein that in humans is encoded by the KIRREL2 gene. Mouse monoclonal antibody raised against a full length recombinant KIRREL2.
Keywords	KIRREL2; kin of IRRE like 2 (Drosophila); NLG1; NEPH3; FILTRIN; kin of IRRE-like protein 2; nephrin-like gene 1; nephrin-like protein 3; kin of irregular chiasm-like protein 2;

GENE INFORMATION

Entrez Gene ID	84063
UniProt ID	Q6UWL6
Pathway	Interactions of the immunoglobulin superfamily (IgSF) member proteins, organism-specific biosystem; Nephrin interactions, organism-specific biosystem