



Magic™ Anti-Pangasius Catfish IgM monoclonal antibody (CABT-B1053)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This monoclonal antibody (Mab) reacts with Anti-Pangasius catfish (Pangasius hypophthalmus) IgM monoclonal antibody immunoglobulin M (IgM) and recognises the heavy chain of the molecule.
Target	Pangasius Catfish (Pangasius hypophthalmus) IgM
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Pangasius Catfish
Conjugate	Unconjugated
Applications	ELISA
Format	Lyophilized
Size	200 µg
Buffer	PBS
Preservative	None
Storage	Store at -20°C prior to reconstitution. For prolonged storage, the Mab solution should be stored at -20°C as working aliquots. Repeated freeze/thawing of the product should be avoided.

BACKGROUND

Keywords

Pangasius Catfish;Pangasius hypophthalmus

GENE INFORMATION

Synonyms

Pangasius Catfish;Pangasius hypophthalmus
