



Mouse anti-Human KLHDC4 monoclonal antibody, clone 5H22 (CABT-B10530)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLHDC4 (AAH01044,2a.a. ~ 463 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5H22
Conjugate	Unconjugated
Applications	WB, IHC, IF, sELISA, ELISA
Sequence Similarities	MGKKKGKKEKKGRGAEKTAAMEKKVSKRSRKEETFLYNELYVYNTRKDTWTKVDIPSPPP RRCahQAVVVPQGGGQLWVFGGEFASPNGEQFYHYKDLWLHLATKTWEQVKSTGGPSGR SGHRMVAWKRQLILFGGFHESTRDYIYYNDVYAFNLDTFTWSKLSPSGTGTPRSGCQMS VTPQGGIVVYGGYSKQRVKKDVKGTRHSDMFLLKPEDGREDKWWTRMNP SGVKPTPRS GFSVAMAPNHQT
Format	Liquid
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords

KLHDC4; kelch domain containing 4; kelch domain-containing protein 4;

GENE INFORMATION

Entrez Gene ID[54758](#)**UniProt ID**[Q8TBB5](#)