



Mouse anti-Human KLHL25 monoclonal antibody, clone 2F0 (CABT-B10531)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLHL25 (AAH28100, 1 a.a. ~ 590 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	2F0
Conjugate	Unconjugated
Applications	WB, ELISA
Sequence Similarities	MSVSVHETRKSRSSSTGSMNVTLFHKASHPDCVLAHLNLTNRKHCMTDVTWAGDRAFPCH RAVLAASSRYFEAMFSHGLRESRDDTVNFDNLHPEVLELLLDFAVSSRIINEENAESL LEAGDMLQFHDVRDAAAEFLEKNLFPSNCLGMMLLSDAHQCRRLYEFSWRMCLVHFETVR QSEDFNSLSKDTLLDLISSDELETEDERVVFEAILQWVKHDLEPRKVHLPPELLRSVRLAL LPSDCLQEAVSSEAL
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Ectoderm-neural cortex protein 2 is a protein that in humans is encoded by the KLHL25 gene. Mouse monoclonal antibody raised against a full length recombinant KLHL25.
Keywords	KLHL25; kelch-like family member 25; ENC2; ENC-2; kelch-like protein 25; kelch-like 25; ectodermal-neural cortex 2; BTB/POZ KELCH domain protein; ectoderm-neural cortex protein 2;

GENE INFORMATION

Entrez Gene ID	64410
UniProt ID	Q9H0H3