



Mouse anti-Human L3MBTL4 monoclonal antibody, clone 2C02D23 (CABT-B10541)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	L3MBTL4 (AAH39316, 1 a.a. ~ 535 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2C02D23
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	MKQPNRKRKLNMDSKERLDQDGRLEQAEKKPKDSTTPLSHVPSAAAQGAWSWEWYLKE QKAVAAPVELFSKDQSFPEHENGFGQIGMRLEGIDPRHPSVFCVLSVAEVCGYRLRLHFDG YLSCYDFWTNAGSPDIHPVGWCEKTKHELHIPKGYRKDKFVWMDYLKACKLQNAPKKLFR NRSPNGPMSKEFQVGMKLEAVDRKNPSLVCVATIADIVEDRLLVHFDNWDDSYDYWCDEVN SPYVQPVGWCQENGR
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	L3MBTL4; l(3)mbt-like 4 (Drosophila); HsT1031; lethal(3)malignant brain tumor-like protein 4; l(3)mbt-like protein 4; H-l(3)mbt-like protein 4;
-----------------	---

GENE INFORMATION

Entrez Gene ID	91133
-----------------------	-----------------------

UniProt ID	Q8NA19
-------------------	------------------------

Function	sequence-specific DNA binding transcription factor activity; zinc ion binding
-----------------	---
