



Mouse anti-Human LAD1 monoclonal antibody, clone 3B0 (CABT-B10542)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LAD1 (NP_005549, 418 a.a. ~ 518 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3B0
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	SESVKSRGLPCTELFVAPVGVASKRHLFEKELAGQSRAEPASSRKENLRSLSGVVTSRLNL WISRTQESGDQDPQEAQKASSATERTQWGQKSDSSLDAEV*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The protein encoded by this gene may be an anchoring filament that is a component of basement membranes. It may contribute to the stability of the association of the epithelial layers
---------------------	---

with the underlying mesenchyme. [provided by RefSeq, Jul 2008]

Keywords

LAD1; ladinin 1; LadA; ladinin-1; lad-1; linear IgA disease antigen homolog; linear IgA bullous dermatosis antigen;

GENE INFORMATION

Entrez Gene ID

[3898](#)

UniProt ID

[O00515](#)

Function

structural molecule activity
