



Mouse anti-Human LBX2 monoclonal antibody, clone 4B6 (CABT-B10550)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LBX2 (NP_001009812, 95 a.a. ~ 195 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	4B6
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	ELERRFVFQKYLAPSERDGLATRLGLANAQVVTWFQNRRAKLKRDVEEMRADVASLRALS PEVLCSLALPEGAPDPGLCLGPAGPDSRPHLSDEEIQVDD*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	LBX2; ladybird homeobox 2; LP3727; transcription factor LBX2; ladybird-like homeobox 2; ladybird-like homeobox 2; ladybird homeobox homolog 2; ladybird homeobox protein homolog 2;
-----------------	---

GENE INFORMATION

Entrez Gene ID [85474](#)

UniProt ID [Q6XYB7](#)

Function sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity
