



Mouse anti-Human LDHAL6B monoclonal antibody, clone 2B4 (CABT-B10554)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LDHAL6B (-, 282 a.a. ~ 382 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	2B4
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	EVTATAYEIIKMKGYTSWAIGLSVADLTESILKNLRRHPVSTIIKGLYGIDEEVFLSIP CILGENGITNLIKIKLTPEEEAHLKKSATLWEIQNKKL*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	LDHAL6B; lactate dehydrogenase A-like 6B; LDHL; LDH6B; LDHAL6; L-lactate dehydrogenase A-like 6B; lactate dehydrogenase A-like 6;
-----------------	---

GENE INFORMATION

Entrez Gene ID	92483
UniProt ID	Q9BYZ2
Pathway	Cysteine and methionine metabolism, organism-specific biosystem; Cysteine and methionine metabolism, conserved biosystem; Glycolysis / Gluconeogenesis, organism-specific biosystem; Glycolysis / Gluconeogenesis, conserved biosystem; Glycolysis and Gluconeogenesis, organism-specific biosystem; Metabolic pathways, organism-specific biosystem
Function	L-lactate dehydrogenase activity; oxidoreductase activity; protein binding