



Mouse anti-Human LEPREL4 monoclonal antibody, clone 2F23 (CABT-B10557)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SC65 (NP_006446, 270 a.a. ~ 380 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	2F23
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	QCKVDCEANLTPNVGGYFVDKFBATMYHYLQFAYYKLNDVRQAARSAASYMLFDPKDSVM QQNLVYYRFHRRARWGLEEEEDFQPREEAMLYHNQTAELELLEFTHMYLQS*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	NO55; SC65; NOL55; synaptonemal complex protein SC65; OTTHUMP00000164758; leprecan-like protein 4; nucleolar autoantigen No55; nucleolar autoantigen (55kD); nucleolar
-----------------	--

GENE INFORMATION

Entrez Gene ID	10609
----------------	-----------------------

UniProt ID	Q92791
------------	------------------------

Function	binding
----------	---------
