



# Rat anti-Mouse Lgi4 monoclonal antibody, clone LU29 (CABT-B10559)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with BSA) corresponding to mouse Lgi4.
Isotype	IgG1
Source/Host	Rat
Species Reactivity	Mouse
Clone	LU29
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS (0.1% sodium azide)
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	LG14; leucine-rich repeat LGI family, member 4; LGIL3; leucine-rich repeat LGI family member 4; LGI1-like protein 3; leucine-rich glioma-inactivated gene 4; leucine-rich glioma-inactivated protein 4;
----------	---

## GENE INFORMATION

**Entrez Gene ID** [243914](#)

---

**Function** molecular\_function;

---