



Mouse anti-Human LMO4 monoclonal antibody, clone 3C7 (CABT-B10576)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LMO4 (AAH17673,2a.a. ~ 165 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	3C7
Conjugate	Unconjugated
Applications	IHC, IF, sELISA, ELISA
Sequence Similarities	MVNP GSSSQPPPV TAGSLSWKRCAGCGGKIADRFLLYAMDSYWHSRCLKSCCQAQLGDI GTSCYTKSGMILCRNDYIRLFGNSGACSACGQSIPASELVMRAQGNVYHLKCFTCSTCRN RLVPGDRFHYINGSLFCEHDRPTALINGHLNSLQSNPLLPDQKVC
Format	Liquid
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This gene encodes a cysteine-rich protein that contains two LIM domains but lacks a DNA-binding homeodomain. The encoded protein may play a role as a transcriptional regulator or as
---------------------	---

an oncogene. [provided by RefSeq, Aug 2008]

Keywords	LMO4; LIM domain only 4; LIM domain transcription factor LMO4; LMO-4; breast tumor autoantigen; LIM domain only protein 4;
-----------------	--

GENE INFORMATION

Entrez Gene ID	8543
-----------------------	----------------------

UniProt ID	P61968
-------------------	------------------------

Pathway	IL6-mediated signaling events, organism-specific biosystem;
----------------	---

Function	enhancer sequence-specific DNA binding; metal ion binding; sequence-specific DNA binding transcription factor activity; transcription factor binding; zinc ion binding;
-----------------	--
