



Mouse anti-Human LMTK3 monoclonal antibody, clone 2F2 (CABT-B10580)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LMTK3 (XP_055866, 1151 a.a. ~ 1250 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2F2
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	VLVNGGLTPPKSEDKVSENGGLRFPNTERPPETGPWRAPGPWEKTPESWGPAPTIGEPA PETSLERAPAPSAVVSSRNGGETAPGPLGPAPKNGTLEPG
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	LMTK3; lemur tyrosine kinase 3; LMR3; TYKLM3; PPP1R101; serine/threonine-protein kinase LMTK3; protein phosphatase 1, regulatory subunit 101;
-----------------	---

GENE INFORMATION

Entrez Gene ID	114783
Function	ATP binding; molecular_function; nucleotide binding; protein serine/threonine kinase activity; protein tyrosine kinase activity