



Mouse anti-Human LNX2 monoclonal antibody, clone 5G8 (CABT-B10582)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LNX2 (NP_699202, 1 a.a. ~ 101 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5G8
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MGTTSDDEMVSVEQTSSSSLNPLCFECGQQHWTRENHLYNYQNEVDDDLVCHICLQPLLQP LDTPCGHTFCYKCLRNFLQEKFDFCPLDRKRLHFKLCKKSS*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	LNX2; ligand of numb-protein X 2; PDZRN1; ligand of Numb protein X 2; numb-binding protein 2; PDZ domain containing ring finger 1; PDZ domain-containing RING finger protein 1;
-----------------	---

GENE INFORMATION

Entrez Gene ID	222484
UniProt ID	Q8N448
Function	PDZ domain binding; metal ion binding; zinc ion binding