



Mouse anti-Human LRIG3 monoclonal antibody, clone 5E6 (CABT-B10588)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LRIG3 (NP_700356, 1020 a.a. ~ 1120 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	5E6
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	LDFSANPEPASVASSNSFMGTFGKALRRPHLDAYSSFGQPSDCQPRAFYLKAHSSPDLDS GSEEDGKERTDFQEENHICTFKQTLENYRTPNPFQSYDLDT*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	LRIG3; leucine-rich repeats and immunoglobulin-like domains 3; LIG3; leucine-rich repeats and immunoglobulin-like domains protein 3; LIG-3;
-----------------	---

GENE INFORMATION

Entrez Gene ID [121227](#)

UniProt ID [Q6UXM1](#)
