



Mouse anti-Human LY86 monoclonal antibody, clone 2I4 (CABT-B10598)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LY86 (AAH38846, 26 a.a. ~ 125 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	2I4
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	AWPTHVVCSDSGLEVLYQSCDPLQDFGFSVEKCSKQLKSNINIRFGIILREDIKELFLDL ALMSQGSSVLNFSYPICEAALPKFSFCGRRKGEQIYYAGP
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	LY86; lymphocyte antigen 86; MD1; MD-1; MMD-1; dJ80N2.1; ly-86; protein MD-1; MD-1, RP105-associated;
-----------------	---

GENE INFORMATION

Entrez Gene ID	9450
UniProt ID	O95711
Pathway	Immune System, organism-specific biosystem; Innate Immunity Signaling, organism-specific biosystem; Toll Like Receptor 4 (TLR4) Cascade, organism-specific biosystem; Toll Receptor Cascades, organism-specific biosystem