



# Magic™ Anti-C. difficile Binary Toxin B Subunit polyclonal antibody, [Biotin] (CABT-B1085)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	C. difficile Binary Toxin B Subunit
Immunogen	C. difficile Binary Toxin, B Subunit
Isotype	IgY
Source/Host	Chicken
Species Reactivity	C. difficile
Purification	Affinity purified
Conjugate	Biotin
Applications	sELISA
Concentration	Store at 2-8°C
Size	100 µg
Buffer	PBS, pH 7.3 with 0.005% thimerosal and 50% glycerol
Preservative	None
Storage	Store at 2-8°C