



Mouse anti-Human UBR3 monoclonal antibody, clone 6B21 (CABT-B11736)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ZNF650 (NP_742067, 1 a.a. ~ 111 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	6B21
Conjugate	Unconjugated
Applications	WB, sELISA, ELISA
Sequence Similarities	MNKRIIEEICRKVTPPVPPKKVTA AEKKTLDKEERRQKARERQQKLLAEFASRQKSFMET AMDVDSPENDIPMEITTAEPQVSEAVYDCVICGQSGPSSDRPTGLVLL*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords UBR3; ubiquitin protein ligase E3 component n-recogin 3 (putative); ZNF650; E3 ubiquitin-protein ligase UBR3; N-recogin-3; zinc finger protein 650; ubiquitin-protein ligase E3-alpha-3;

ubiquitin-protein ligase E3-alpha-III;

GENE INFORMATION

Entrez Gene ID [130507](#)

UniProt ID [Q6ZT12](#)

Function ligase activity; metal ion binding; ubiquitin-protein ligase activity; zinc ion binding
