



Anti-Casp2 (N-terminal) polyclonal antibody (CABT-B438)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Casp2
Immunogen	Synthetic peptide derived from the N-terminal region of zebrafish Caspase-2 protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Zebrafish
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µl
Buffer	MOPS buffer (pH 7.5) containing 0.1% BSA, 0.05% NaN ₃ and 50% glycerol
Preservative	0.05% Sodium Azide
Storage	Store at -20°C for up to 24 months upon receiving the product.

BACKGROUND

Keywords	CASP2;caspase 2, apoptosis-related cysteine peptidase;ICH1;NEDD2;CASP-2;NEDD-2;PPP1R57;caspase-2;protease ICH-1;protein phosphatase 1, regulatory subunit 57;neural precursor cell expressed developmentally down-regulated protein 2
-----------------	---

GENE INFORMATION

Synonyms

casp2; caspase 2, apoptosis-related cysteine protease
