



Anti-3-Mar polyclonal antibody (CABT-B704)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	3-Mar
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Dog, Pig, Rat, Zebrafish, Guinea pig, Horse, Rabbit, Bovine
Conjugate	Unconjugated
Applications	WB
Format	Lyophilized powder
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration)
Size	50 µg
Buffer	PBS buffer with 2% sucrose
Preservative	None
Storage	For longer periods of storage, store at -20°C. Avoid repeat freeze-thaw cycles.

BACKGROUND

Keywords MARCH3;membrane-associated ring finger (C3HC4) 3;E3 ubiquitin-protein ligase
MARCH3;FLJ20096;MARCH-III;Membrane associated ring finger (C3HC4) 3;Membrane-associated RING finger protein 3;Membrane-associated RING-CH protein III;MGC48332

GENE INFORMATION

Synonyms

March3; 6330411I15Rik; MARCH-III; A530081L18Rik; membrane-associated ring finger (C3HC4) 3; AI843639
