



## Anti-Mpx polyclonal antibody (CABT-B738)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Target</b>	Mpx
<b>Immunogen</b>	Recombinant protein encompassing a sequence within the center region of zebrafish Mpx. The exact sequence is proprietary.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Zebrafish
<b>Purification</b>	Purified by antigen-affinity chromatography.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-Wmt
<b>Format</b>	Liquid
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration)
<b>Size</b>	50 µl
<b>Buffer</b>	1XPBS, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
<b>Preservative</b>	None
<b>Storage</b>	Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

### BACKGROUND

**Keywords**

MPX;myeloid-specific  
peroxidase;ENSDARG00000019521;mpo;fj80f04;wu:fj80f04;chromosome:  
10;eosinophil peroxidase;eosinophil peroxidase;myeloperoxidase;EC  
1.11.1.7;

---

## GENE INFORMATION

**Synonyms**

mpx; fj80f04; myeloid-specific peroxidase; mpo; wu:fj80f04; myeloperoxidase

---