



Magic™ Anti-NSE monoclonal antibody, clone N23720OT2 (CABT-B8428)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human NSE
Immunogen	NSE antibody was raised in mouse using purified NSE from human brain as the immunogen.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	N23720OT2
Purification	Protein G affinity chromatography
Conjugate	Unconjugated
Applications	ELISA (cap), IHC, IA, WB
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4, 0.1 % sodium azide (NaN ₃)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C