



Anti-CA 24-2 monoclonal antibody, clone N73973 (CABT-B8458)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Human colorectal adenocarcinoma cell-line COLO-205
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	N73973
Purification	Protein A affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Buffer	0.01M PBS buffer,pH 7.4
Preservative	None
Storage	Store at 4°C for short term storage, store at -20°C or -80°C for long term storage

BACKGROUND

Introduction	The CA-242 is a sialylated carbohydrate antigen present on mucinous type of glycoproteins in carcinomas of many organs. The CA-242 antigen is shedded from the tumor and the CA-242
--------------	---

can be detected in serum from patients with carcinomas.
