



# Magic™ Anti-Malaria (pan-pLDH) Monoclonal Antibody, clone IN361 (CABT-B8661)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Malaria pLDH. Detection capability of $\geq 50$ parasites/ml
<b>Immunogen</b>	Recombinant Plasmodium LDH
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Plasmodium
<b>Purification</b>	Ion Exchange Purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA(cap), LFIA
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.