



Magic™ Anti-Malaria (pan-pLDH) Monoclonal Antibody, clone IN360 (CABT-B8662)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Malaria pLDH. Detection capability of ≥ 50 parasites/ml
Immunogen	Recombinant Plasmodium LDH
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Plasmodium
Purification	Ion Exchange Purified
Conjugate	Unconjugated
Applications	ELISA(det), LFIA
Format	Liquid
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Preservative	0.1% Sodium Azide
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.