



Anti-Malaria Aldolase polyclonal Antibody (CABT-B8664)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	highly purified recombinant plasmodium aldolase antigen
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Plasmodium
Conjugate	Unconjugated
Applications	ELISA, WB, IP
Format	Liquid
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C