



Anti-NDV HN protein polyclonal antibody [HRP] (CABT-B8668)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide derived from NDV HN protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	NDV
Purification	Affinity purified by Protein A
Conjugate	HRP
Applications	ELISA, WB, IHC-P, IHC-Fr
Format	Lyophilized or Liquid
Size	100 µl
Buffer	0.01M TBS(pH7.4) with 1% BSA and 50% Glycerol.
Preservative	None
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of anti