



Anti-SEZ6L polyclonal antibody (CABT-B8670)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide derived from human SEZ6L
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified by Protein A
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P, IHC-Fr, IF
Format	Liquid
Concentration	1.0 mg/ml
Buffer	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Preservative	None
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C.

BACKGROUND

Introduction	SEZ6L is a candidate tumor suppressor gene. This protein may contribute to specialized endoplasmic reticulum functions in neurons.
--------------	--