



Anti-NDV HN protein polyclonal antibody [Alexa Fluor® 488] (CABT-B8672)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide derived from NDV HN protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	NDV
Purification	Affinity purified by Protein A
Conjugate	Alexa Fluor 488
Applications	IF
Format	Lyophilized or Liquid
Buffer	1% BSA, 0.01M PBS, pH 7.4.
Preservative	15mM Sodium Azide
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of anti