



Anti-Rubella virus Polyclonal antibody [FITC] (CABT-B8675)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with Rubella virus (Purified virions)
Immunogen	strain HPV77
Isotype	IgG
Source/Host	Goat
Species Reactivity	Rubella virus
Purification	Purified by protein A chromatography or sequential differential precipitations
Conjugate	FITC
Applications	ICC/IF, HI
Format	Liquid
Buffer	Phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative
Preservative	0.1% Sodium Azide
Storage	Store at -20°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C. at -20°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C.