



Anti-CMV IE72 monoclonal antibody, clone EG23 (CABT-B8688)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cytomegalovirus Immediate Early Protein
Isotype	IgG3
Source/Host	Mouse
Species Reactivity	Cytomegalovirus
Clone	EG23
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IF, ELISA, WB
Format	Liquid
Concentration	1 mg/mL
Buffer	Phosphate Buffered Saline pH 7.4
Preservative	None
Storage	Upon receipt, store at -20°C or below. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antibody be stored in use size aliquots and thawed just prior to use. Storage at 2 - 8°C for up to two weeks is acceptable fo