



Magic™ Anti-E. coli Enterotoxin STp monoclonal antibody, clone C2037N (CABT-B8713)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to E. coli Heat-Stable Enterotoxin (STp).
Immunogen	STp peptide (18 a.a. length).
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	E. coli
Clone	C2037N
Purification	Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA (cap)
Format	Liquid
Concentration	Lot specific
Buffer	PBS: 3.0 mM Monopotassium Phosphate (KH ₂ PO ₄), 8.0 mM Disodium Phosphate (Na ₂ HPO ₄), 1.5 mM Potassium Chloride (KCl), 140 mM Sodium Chloride (NaCl), pH 7.4.
Preservative	0.05% Sodium Azide
Storage	Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple

freeze/thaw cycles.
