



# Anti-S. pneumoniae CWPS monoclonal antibody (CABT-B8727)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Streptococcus pneumoniae common cell wall polysaccharide (CWPS).
<b>Immunogen</b>	Purified Streptococcus pneumoniae CWPS.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	S. pneumoniae
<b>Purification</b>	Protein A Chromatograph
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA (det), IF, LFIA
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	0.01 M Phosphate Buffered Saline, pH 7.2
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 2-8°C.