



Anti-EBV EBNA-1 polyclonal antibody (CABT-B8734)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	EBNA-1
Immunogen	Purified recombinant EBNA-1
Isotype	IgG
Source/Host	Goat
Species Reactivity	EBV
Conjugate	Unconjugated
Applications	ELISA, IF, WB
Format	Liquid
Size	500 µl
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.