



## Anti-Coxsackievirus A6 VP1 polyclonal antibody (CABT-B8738)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein corresponding to Coxsackievirus A6 virus VP1 protein. The exact sequence is proprietary.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Coxsackievirus A6
<b>Purification</b>	Purified by antigen-affinity chromatography.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration)
<b>Size</b>	100 µl
<b>Buffer</b>	1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.
<b>Preservative</b>	None
<b>Storage</b>	Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.