



Anti-Neuropeptide Y polyclonal antibody (CABT-B8739)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic Neuropeptide Y (porcine) conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Monkey, Pig, Rat, Cat, Rabbit, Bovine
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Dot, ELISA, IHC-Fr, RIA
Format	Liquid
Size	100 µl, 200 µl
Buffer	Whole antiserum containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Keep as concentrated solution. Store at 4°C short term. For extended storage aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

BACKGROUND

Introduction	NPY, a 36 amino acid peptide amide isolated from porcine brain, is a major regulatory neuropeptide in the mammalian central and peripheral nervous systems. NPY belongs to the pancreatic polypeptide family of peptides which are characterized by a common tertiary
---------------------	---

structure.
