



Anti-SARS-CoV Nucleoprotein monoclonal antibody, clone CEJ297 (CABT-B8748)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific for the nucleoprotein of the SARS virus. Does not react with human coronaviruses 229E and OC43 as well as FIP-1, FIP 2, canine coronavirus, TGEV and mouse hepatitis virus.
Immunogen	Recombinant SARS Nucleoprotein
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	SARS
Clone	CEJ297
Purification	Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA, WB, IF
Format	Liquid
Size	1 mg
Buffer	0.01 M PBS, pH 7.2 This product contains no stabilizing proteins.
Preservative	0.1% Sodium Azide
Storage	Upon receipt, aliquot and store at –20°C. Avoid multiple freeze/thaw cycles.