



Magic™ Anti-SAA (aa 23-29) monoclonal antibody (CABT-B8775)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes native and recombinant Human SAA. Cross reacts with canine, equine and feline SAA.
Immunogen	Synthetic peptides representing the region aa 23-29 of human SAA
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human, Canine, Equine and Feline
Purification	Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA(cap), WB
Format	Liquid
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C.