



Magic™ Anti-SAA monoclonal antibody (CABT-B8776)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes native and recombinant Human SAA. Cross reacts with canine, equine and feline.
Immunogen	Canine SAA
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human, Canine, Equine and Feline
Purification	Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA(det), WB
Format	Liquid
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C.