



# Anti-HPV6 L1 polyclonal antibody (CABT-B8785)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                      |
|--------------------|------------------------------------------------------|
| Isotype            | IgG                                                  |
| Source/Host        | Rabbit                                               |
| Species Reactivity | HPV6                                                 |
| Purification       | Protein G purified                                   |
| Conjugate          | Unconjugated                                         |
| Applications       | WB, ELISA                                            |
| Format             | Liquid                                               |
| Concentration      | 1mg/1ml                                              |
| Size               | 100 µl                                               |
| Buffer             | Tris-Gly                                             |
| Preservative       | None                                                 |
| Storage            | Store at -20°C. Avoid repeated freezing and thawing. |