



Anti-HPV11 L1 polyclonal antibody (CABT-B8787)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|--|
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | HPV11 |
| Purification | Protein G purified |
| Conjugate | Unconjugated |
| Applications | WB, ELISA |
| Format | Liquid |
| Concentration | 1mg/1ml |
| Size | 100 µg |
| Buffer | Tris-Gly |
| Preservative | None |
| Storage | Store at -20°C. Avoid repeated freezing and thawing. |