



Anti-HPV45 L1 polyclonal antibody (CABT-B8797)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HPV45
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	1mg/1ml
Size	100 µl
Buffer	Tris-Gly
Preservative	None
Storage	Store at -20°C. Avoid repeated freezing and thawing.