



## Anti-HPV52 L1 polyclonal antibody (CABT-B8799)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	HPV52
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Format</b>	Liquid
<b>Concentration</b>	1mg/1ml
<b>Size</b>	100 µg
<b>Buffer</b>	Tris-Gly
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. Avoid repeated freezing and thawing.