



# Anti-TMEM26 polyclonal antibody (CABT-B8820)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A portion of amino acids 300-350 of human Tmem26 was used as the immunogen.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, FC
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	PBS
<b>Preservative</b>	0.05% Sodium Azide
<b>Storage</b>	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

## BACKGROUND

<b>Introduction</b>	Transmembrane protein 26, or Tmem26, is a protein which spans the entire biological membrane and thus maintains three dimensional structures. Transmembrane proteins aggregate and precipitate in water. They require detergents or nonpolar solvents for extraction; although some of them (beta-barrels) can be also extracted using denaturing agents.
---------------------	---

# GENE INFORMATION

Entrez Gene ID [219623](#)

---

UniProt ID [Q6ZUK4](#)

---