



Anti-MBTPS1 polyclonal antibody (CABT-B8877)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein of human MBTPS1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	8720
UniProt ID	Q14703