



Anti-Triiodothyronine polyclonal antibody (CABT-B8931)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	T3-BgG (carboxyl derivative)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Protein A Purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Preservative	0.05% Sodium Azide
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

BACKGROUND

Introduction	Triiodothyronine, C ₁₅ H ₁₂ I ₃ NO ₄ , also known as T ₃ , is a thyroid hormone. It affects almost every physiological process in the body, including growth and development, metabolism, body temperature, and heart rate.
--------------	--